

HOLLYWOOD SMILES DENTAL LABORATORY LTD.

Therm-O-Flex™ A new night guard has arrived.

Therm-O-Flex™ a thermo activated material, is a unique product compared to the typical hard acrylic material. Night guards are typically made from a hard acrylic material. A material that could easily be broken and usually causes discomfort for the patient. Now, with **Therm-O-Flex™** night guards from Hollywood Smiles Dental Laboratory, a new and truly unique product has arrived.

Unique Therm-O-Flex™ Material

- resistant to breakage
- reduces chair time
- requires minimal adjustments, if any
- softens slightly in the mouth for added comfort

Accurate Fit

Softens slightly in warm water, or in the oral cavity, thus conforming more readily to tooth contours.

Memory

When distorted, it returns to the processed state, thus maintaining fit integrity.

Retention

The thermo elasticity allows it to gain retention by engaging the lingual and buccal interproximal spaces. It also maintains retention from the lingual and buccal undercuts under height of contour of teeth.

Prescription Requirements

Upper and lower models and a bite registration is required.

